

Antibodies Product Datasheet

Product Name:	ERCC6L2 Antibody
Catalog Number:	C50245
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ERCC6L2
Concentration:	1 mg/ml
Gene Name:	RAD26L C9orf102
Protein Name:	Putative DNA repair and recombination protein RAD26-like (EC 3.6.4.-)
Other Name:	
Human Swiss Prot No.:	Q5T890
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ERCC6L2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year