

Antibodies Product Datasheet

Product Name:	Cytokeratin 8 Monoclonal Antibody
Catalog Number:	N0189
Reactivity:	Human
Applications:	WB,IHC-p,IF/ICCELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human Cytokeratin 8 (aa391-483) expressed in E. Coli.
Concentration:	
Gene Name:	KRT8
Protein Name:	Keratin type II cytoskeletal 8
Other Name:	KRT8; CYK8; Keratin; type II cytoskeletal 8; Cytokeratin-8; CK-8; Keratin-8; K8; Type-II keratin Kb8
Human Swiss Prot No.:	P05787
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Cytokeratin 8 Monoclonal Antibody detects endogenous levels of Cytokeratin 8 protein.
Storage and Stability:	-20°C/1 year