

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | ACE2 Antibody |
| Catalog Number: | C42051 |
| Reactivity: | Human, Mouse, Rat |
| Applications: | IHC-p, WB |
| Dilutions: | IHC-p1:50-200 ,WB 1:1000-2000 |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Immunogen: | Synthesized peptide derived from human ACE2 |
| Concentration: | 1 mg/ml |
| Gene Name: | ACE2 UNQ868/PRO1885 |
| Protein Name: | ACE2 |
| Other Name: | Angiotensin-converting enzyme 2 (EC 3.4.17.23) (ACE-related carboxypeptidase) (Angiotensin-converting enzyme homolog) (ACEH) (Metalloprotease MPROT15) [Cleaved into: Processed angiotensin-converting enzyme 2]; YM33046 |
| Human Swiss Prot No.: | Q9BYF1 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.06% sodium azide. |
| Specificity: | This antibody detects endogenous levels of ACE2 at Human, Mouse,Rat |
| Storage and Stability: | -20°C/1 year |