

Antibodies Product Datasheet

Product Name:	Smad1/5/9 (Phospho-Ser463+Ser465) Antibody
Catalog Number:	A8417
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Smad1/5/9 (Phospho-Ser463+Ser465)
Concentration:	1 mg/ml
Gene Name:	SMAD1 BSP1 MADH1 MADR1
Protein Name:	Smad1/5/9 (Phospho-Ser463+Ser465)
Other Name:	Mothers against decapentaplegic homolog 1 (MAD homolog 1) (Mothers against DPP homolog 1) (JV4-1) (Mad-related protein 1) (SMAD family member 1) (SMAD 1) (Smad1) (hSMAD1) (Transforming growth factor-beta-signaling protein 1) (BSP-1)
Human Swiss Prot No.:	Q15797
Formulation:	Liquid in PBS containing 50% glycerol, and 0.188% sodium azide.
Specificity:	This antibody detects endogenous levels of Smad1/5/9 (Phospho-Ser463+Ser465) at Human, Mouse, Rat
Storage and Stability:	-20°C/1 year