

## Antibodies Product Datasheet

<b>Product Name:</b>	FAK Antibody
<b>Catalog Number:</b>	C42198
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	IHC-p, WB
<b>Dilutions:</b>	IHC-p1:50-200 ,WB 1:1000-2000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human FAK
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	PTK2 FAK FAK1
<b>Protein Name:</b>	FAK
<b>Other Name:</b>	Focal adhesion kinase 1 (FADK 1) (EC 2.7.10.2) (Focal adhesion kinase-related nonkinase) (FRNK) (Protein phosphatase 1 regulatory subunit 71) (PPP1R71) (Protein-tyrosine kinase 2) (p125FAK) (pp125FAK); YM3814
<b>Human Swiss Prot No.:</b>	Q05397
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.86% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of FAK at Human, Mouse, Rat
<b>Storage and Stability:</b>	-20°C/1 year