

Antibodies Product Datasheet

Product Name:	Cytokeratin 18 Monoclonal Antibody(8F2)
Catalog Number:	N1362
Reactivity:	Human,Mouse,Rat
Applications:	WB
Dilutions:	WB: 1:1000-3000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of Keratin type I cytoskeletal 18
Concentration:	
Gene Name:	KRT18
Protein Name:	Keratin type I cytoskeletal 18
Other Name:	KRT18; CYK18; PIG46; Keratin, type I cytoskeletal 18; Cell proliferation-inducing gene 46 protein; Cytokeratin-18; CK-18; Keratin-18; K18
Human Swiss Prot No.:	P05783
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Cytokeratin 18 protein.
Storage and Stability:	-20°C/1 year