

Antibodies Product Datasheet

Product Name:	I κ B- α Antibody
Catalog Number:	C10481
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, WB, IHC-p, ELISA
Dilutions:	IF: 1:50-200 WB 1:500-2000, ELISA 1:10000-20000 IHC 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human I κ B- α . AA range:1-50
Concentration:	1 mg/ml
Gene Name:	NFKBIA
Protein Name:	NF-kappa-B inhibitor alpha
Other Name:	NFKBIA; IKBA; MAD3; NFKBI; NF-kappa-B inhibitor alpha; I-kappa-B-alpha; I κ B-alpha; I κ B α ; Major histocompatibility complex enhancer-binding protein MAD3
Human Swiss Prot No.:	P25963
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	I κ B- α Polyclonal Antibody detects endogenous levels of I κ B- α protein.
Storage and Stability:	-20°C/1 year