

## Antibodies Product Datasheet

<b>Product Name:</b>	PEG10 mouse mAb
<b>Catalog Number:</b>	N0919
<b>Reactivity:</b>	Human
<b>Applications:</b>	WB,IP
<b>Dilutions:</b>	wb 1:1000
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Purified recombinant human PEG10 protein fragments expressed in E.coli.
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	peg10
<b>Protein Name:</b>	
<b>Other Name:</b>	AA407948;Edr;Embryonal carcinoma differentiation regulated;Embryonal carcinoma differentiation-regulated protein;HB1;HB1;KIAA1051;Mammalian retrotransposon-derived protein 2;Mar2;Mart2;MEF3 like 1;MEF3-like protein 1;MEF3L;MEF3L1;MyEF3;Myelin expression
<b>Human Swiss Prot No.:</b>	Q86TG7
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of PEG10 and does not cross-react with related proteins.
<b>Storage and Stability:</b>	-20°C/1 year