

Antibodies Product Datasheet

Product Name:	REQU Antibody
Catalog Number:	R12-3421
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	REQU Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human REQU.
Concentration:	1 mg/ml
Gene Name:	DPF2
Other Name:	Zinc finger protein ubi-d4; Apoptosis response zinc finger protein; BRG1-associated factor 45D; BAF45D; D4 zinc and double PHD fingers family 2; Protein requiem; DPF2; BAF45D; REQ; UBID4
Human Swiss Prot No.:	Q92785
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	REQU Antibody antibody detects endogenous levels of REQU.
Storage and Stability:	Stable at -20°C for at least 1 year.