

Antibodies Product Datasheet

Product Name:	p47-phox Antibody
Catalog Number:	B1160
Reactivity:	Human, Monkey
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human p47 phox. AA range:281-330
Concentration:	1 mg/ml
Gene Name:	NCF1
Protein Name:	Neutrophil cytosol factor 1
Other Name:	NCF1; NOXO2; SH3PXD1A; Neutrophil cytosol factor 1; NCF-1; 47 kDa autosomal chronic granulomatous disease protein; 47 kDa neutrophil oxidase factor; NCF-47K; Neutrophil NADPH oxidase factor 1; Nox organizer 2; Nox-organizing protein 2; SH3
Human Swiss Prot No.:	P14598
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	p47-phox Polyclonal Antibody detects endogenous levels of p47-phox protein.
Storage and Stability:	-20°C/1 year