

Antibodies Product Datasheet

Product Name:	IRF-3 (Phospho-Ser386) Antibody
Catalog Number:	A8203
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IRF-3 around the phosphorylation site of Ser386. AA range:352-401
Concentration:	1 mg/ml
Gene Name:	IRF3
Protein Name:	Interferon regulatory factor 3
Other Name:	IRF3; Interferon regulatory factor 3; IRF-3
Human Swiss Prot No.:	Q14653
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-IRF-3 (S386) Polyclonal Antibody detects endogenous levels of IRF-3 protein only when phosphorylated at S386.
Storage and Stability:	-20°C/1 year