

Antibodies Product Datasheet

Product Name:	TRAP220 Antibody
Catalog Number:	B0556
Reactivity:	Human, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PPAR-BP. AA range:1423-1472
Concentration:	1 mg/ml
Gene Name:	MED1
Protein Name:	Mediator of RNA polymerase II transcription subunit 1
Other Name:	MED1; ARC205; CRSP1; CRSP200; DRIP205; DRIP230; PBP; PPARBP; PPARGBP; RB18A; TRAP220; TRIP2; Mediator of RNA polymerase II transcription subunit 1; Activator-recruited cofactor 205 kDa component; ARC205; Mediator complex subunit 1; Peroxiso
Human Swiss Prot No.:	Q15648
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TRAP220 Polyclonal Antibody detects endogenous levels of TRAP220 protein.
Storage and Stability:	-20°C/1 year