

Antibodies Product Datasheet

Product Name:	GGT1(heavy chain, Cleaved-Gly380) Antibody
Catalog Number:	L0452
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human GGT1 (heavy chain, Cleaved-Gly380)
Concentration:	1 mg/ml
Gene Name:	GGT1 GGT
Protein Name:	GGT1 (heavy chain, Cleaved-Gly380)
Other Name:	Gamma-glutamyltranspeptidase 1 (GGT 1) (EC 2.3.2.2) (Gamma-glutamyltransferase 1) (Glutathione hydrolase 1) (EC 3.4.19.13) (Leukotriene-C4 hydrolase) (EC 3.4.19.14) (CD antigen CD224) [Cleaved into: Gamma-glutamyltranspeptidase 1 heavy chain; Gamma-glutam
Human Swiss Prot No.:	P19440
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human GGT1 (heavy chain, Cleaved-Gly380, protein was cleaved amino acid sequence between 380-381)
Storage and Stability:	-20°C/1 year