

Antibodies Product Datasheet

Product Name:	Anti-Cytokeratin 15 [ABT-CK15]
Catalog Number:	V0047
Clone ID:	ABT-CK15
Reactivity:	Human
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	P19012
Gene Name:	KRT15
Other Name:	Keratin, type I cytoskeletal 15; Cytokeratin-15; CK-15; K15
Isotype	Mouse IgG
Localization:	N/A
Positive Tissue Control:	N/A
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	Anti-Cytokeratin 15 [ABT-CK15] detects endogenous levels of total Cytokeratin 15 protein.
Storage and Stability:	-20°C/1 year