

Antibodies Product Datasheet

Product Name:	AF-10 Antibody
Catalog Number:	C10737
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AF10. AA range:211-260
Concentration:	1 mg/ml
Gene Name:	MLLT10
Protein Name:	Protein AF-10
Other Name:	MLLT10; AF10; Protein AF-10; ALL1-fused gene from chromosome 10 protein
Human Swiss Prot No.:	P55197
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AF-10 Polyclonal Antibody detects endogenous levels of AF-10 protein.
Storage and Stability:	-20°C/1 year