

Antibodies Product Datasheet

Product Name:	TTC7B Antibody
Catalog Number:	C50322
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human TTC7B
Concentration:	1 mg/ml
Gene Name:	TTC7B TTC7L1
Protein Name:	Tetratricopeptide repeat protein 7B (TPR repeat protein 7B) (Tetratricopeptide repeat protein 7-like-1) (TPR repeat protein 7-like-1)
Other Name:	
Human Swiss Prot No.:	Q86TV6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of TTC7B at Human
Storage and Stability:	Store at -20° C/1 year