

## Antibodies Product Datasheet

<b>Product Name:</b>	Integrin $\beta$ 1 Antibody
<b>Catalog Number:</b>	B0445
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, IHC-p, IF(paraffin section), ELISA
<b>Dilutions:</b>	IHC-p: 100-300. Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human Integrin beta1. AA range:749-798
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	ITGB1
<b>Protein Name:</b>	Integrin beta-1
<b>Other Name:</b>	ITGB1; FN1R; MDF2; MSK12; Integrin beta-1; Fibronectin receptor subunit beta; VLA-4 subunit beta; CD antigen CD29
<b>Human Swiss Prot No.:</b>	P05556
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Integrin $\beta$ 1 Polyclonal Antibody detects endogenous levels of Integrin $\beta$ 1 protein.
<b>Storage and Stability:</b>	-20°C/1 year