

Antibodies Product Datasheet

Product Name:	NKX3-2 Antibody
Catalog Number:	C50045
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human NKX3-2
Concentration:	1 mg/ml
Gene Name:	NKX3-2 BAPX1 NKX3B
Protein Name:	Homeobox protein Nkx-3.2 (Bagpipe homeobox protein homolog 1) (Homeobox protein NK-3 homolog B)
Other Name:	
Human Swiss Prot No.:	P78367
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of NKX3-2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year