

Antibodies Product Datasheet

Product Name:	Acetyl-Histone H3 (K9) Monoclonal Antibody(2E7)
Catalog Number:	N1625
Reactivity:	Human,Rat,Mouse
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of Acetyl-Histone H3 (K9)
Concentration:	1 mg/ml
Gene Name:	HIST1H3A
Protein Name:	HIST1H3A
Other Name:	HIST1H3A;H3k9AC
Human Swiss Prot No.:	P68431
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Acetyl-Histone H3 (K9) Monoclonal Antibody(2E7) detects endogenous levels of HIST1H3A
Storage and Stability:	-20°C/1 year