

Antibodies Product Datasheet

Product Name:	CYP26A1 Antibody
Catalog Number:	C12254
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Cytochrome P450 26A1. AA range:251-300
Concentration:	1 mg/ml
Gene Name:	CYP26A1
Protein Name:	Cytochrome P450 26A1
Other Name:	CYP26A1; CYP26; P450RAI1; Cytochrome P450 26A1; Cytochrome P450 retinoic acid-inactivating 1; Cytochrome P450RAI; hP450RAI; Retinoic acid 4-hydroxylase; Retinoic acid-metabolizing cytochrome
Human Swiss Prot No.:	O43174
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CYP26A1 Polyclonal Antibody detects endogenous levels of CYP26A1 protein.
Storage and Stability:	-20°C/1 year