

Antibodies Product Datasheet

Product Name:	SNRPN Antibody
Catalog Number:	C30295
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human SNRPN. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	SNRPN
Protein Name:	Small nuclear ribonucleoprotein-associated protein N
Other Name:	SNRPN; HCERN3; SMN; Small nuclear ribonucleoprotein-associated protein N; snRNP-N; Sm protein D; Sm-D; Sm protein N; Sm-N; SmN; Tissue-specific-splicing protein
Human Swiss Prot No.:	P63162
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SNRPN Polyclonal Antibody detects endogenous levels of SNRPN protein.
Storage and Stability:	-20°C/1 year