

Antibodies Product Datasheet

Product Name:	Hsp27 mouse mAb
Catalog Number:	N0913
Reactivity:	Human,Monkey
Applications:	WB,ICC
Dilutions:	wb 1:1000 icc 1:300
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human Hsp27 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	hspb1
Protein Name:	
Other Name:	Heat shock 27kDa protein;28 kDa heat shock protein;CMT2F;DKFZp586P1322;Estrogen regulated 24 kDa protein;Estrogen-regulated 24 kDa protein;Heat shock 25kDa protein 1;Heat shock 27 kDa protein;Heat shock 27kD protein 1;Heat shock 27kDa protein 1;Heat shock
Human Swiss Prot No.:	P04792
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Hsp27 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year