
Antibodies Product Datasheet

Product Name:	NF-kappaB p65 Antibody
Catalog Number:	R12-3091
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC, ELISA
Dilutions:	NF-kappaB p65 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NF-kappaB p65.
Concentration:	1 mg/ml
Gene Name:	RELA
Other Name:	NFKB3; nuclear factor NF-kappa-B p65 subunit; p65; RELA; TF65; transcription factor p65
Human Swiss Prot No.:	Q04206
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	NF-kappaB p65 Antibody antibody detects endogenous levels of NF-kappaB p65.
Storage and Stability:	Stable at -20°C for at least 1 year.