

Antibodies Product Datasheet

Product Name:	Anti-HumanPlacental lactogen (hPL) + Growth Hormone (GH) (ABT540) Mouse Monoclonal Antibody
Catalog Number:	V0079
Clone ID:	ABT540
Reactivity:	Human
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	CSH1 CSH2 GH1
Gene Name:	P01241 P0DML2 P0DML3
Other Name:	Lactogen, Placental lactogen, PL, hPL, gH, GH-N, GH1, GHB5, GHN, Growth hormone 1, Growth hormone, Growth hormone B5, Growth hormone, normal, Growth hormone, pituitary, HG1, hGH-N, IGHD1B, Pituitary growth hormone, RNGHGP, SOMA_HUMAN, Somatotropin
Isotype	IgG1, Kappa
Localization:	Pituitary
Positive Tissue Control:	Cytoplasmic
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can recognize human Placental lactogen, and has cross-reactiion with GH protein.
Storage and Stability:	-20°C/1 year