

## Antibodies Product Datasheet

<b>Product Name:</b>	Anti-Human Melan-A (ABT-Mart1) Mouse Monoclonal Antibody
<b>Catalog Number:</b>	V0161
<b>Clone ID:</b>	ABT-Mart1
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	IHC
<b>Application Dilutions</b>	IHC 1:100-200
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>UniProt:</b>	MLANA
<b>Gene Name:</b>	Q16655
<b>Other Name:</b>	Antigen LB39 AA, Antigen LB39-AA, Antigen SK29 AA, Antigen SK29-AA, MAR1_HUMAN, MART 1, MART-1, MART1, Melan A, Melan A protein, Melanoma antigen recognized by T cells 1, Melanoma antigen recognized by T-cells 1, MLAN A, MLANA, OTTHUMP00000021036, OTTHUMP00000021037, OTTHUMP00000021038, Protein Melan-A
<b>Isotype</b>	IgG2b, Kappa
<b>Localization:</b>	Skin/ Melanoma
<b>Positive Tissue Control:</b>	Cytoplasmic
<b>Immunogen</b>	Synthetic peptide
<b>Purification</b>	Immunogen affinity purified
<b>Formulation:</b>	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
<b>Specificity:</b>	The antibody can specifically recognize human Melan-A protein.
<b>Storage and Stability:</b>	-20°C/1 year