

## Antibodies Product Datasheet

<b>Product Name:</b>	IFIT3 Antibody
<b>Catalog Number:</b>	C44200
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	WB 1:500-2000 ELISA 1:5000-20000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human protein . at AA range: 330-410
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	IFIT3 CIG-49 IFI60 IFIT4 ISG60
<b>Protein Name:</b>	Interferon-induced protein with tetratricopeptide repeats 3 (IFIT-3) (CIG49) (ISG-60) (Interferon-induced 60 kDa protein) (IFI-60K) (Interferon-induced protein with tetratricopeptide repeats 4) (IFIT-4) (Retinoic acid-induced gene G protein) (P60) (RIG-G)
<b>Other Name:</b>	
<b>Human Swiss Prot No.:</b>	O14879
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
<b>Specificity:</b>	IFIT3 Polyclonal Antibody detects endogenous levels of protein.
<b>Storage and Stability:</b>	-20°C/1 year