

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | 53BP1 (Phospho-Ser1778) Antibody |
| Catalog Number: | A9771 |
| Reactivity: | Human, Rat, Mouse |
| Applications: | WB |
| Dilutions: | WB 1:1000-2000 |
| Modification: | Phospho Specific |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Immunogen: | Synthesized phospho peptide around human 53BP1 (Ser1778) |
| Concentration: | 1 mg/ml |
| Gene Name: | TP53BP1 |
| Protein Name: | 53BP1 (Ser1778) |
| Other Name: | Tumor suppressor p53-binding protein 1 (53BP1) (p53-binding protein 1) (p53BP1) |
| Human Swiss Prot No.: | Q12888 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | This antibody detects endogenous levels of Human 53BP1 (phospho-Ser1778) |
| Storage and Stability: | -20°C/1 year |