

Antibodies Product Datasheet

Product Name:	Anti-Human CD7 (ABT035) Mouse Monoclonal Antibody
Catalog Number:	V0252
Clone ID:	ABT035
Reactivity:	Human
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	CD7
Gene Name:	P09564
Other Name:	CD7, CD7 antigen (p41), CD7 antigen, CD7 molecule, CD7_HUMAN, GP40, LEU 9, LEU9, p41 protein, T cell antigen CD7, T cell leukemia antigen, T cell surface antigen Leu 9, T-cell antigen CD7, T-cell leukemia antigen, T-cell surface antigen Leu-9, Tp 40, Tp40, TP41
Isotype	IgG2a, Kappa
Localization:	Tonsil/ Appendix
Positive Tissue Control:	Membranous
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human CD7 protein .In western blotting of Jurkat cell lysate, there shows a 37 kDa band detected by the antibody.
Storage and Stability:	-20°C/1 year