

Antibodies Product Datasheet

Product Name:	ERK1/2 mouse Monoclonal Antibody(1H4)
Catalog Number:	N1659
Reactivity:	Human,Rat,Mouse
Applications:	WB,IHC-p,IF(paraffin section)
Dilutions:	WB 1:1000-2000, IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of ERK1/2 at AA range of 140-220
Concentration:	1 mg/ml
Gene Name:	MAPK1/MAPK3
Protein Name:	MAPK1/MAPK3
Other Name:	MAPK1/MAPK3
Human Swiss Prot No.:	P27361/P28482
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ERK1/2 protein detects endogenous levels of MAPK1/MAPK3
Storage and Stability:	-20°C/1 year