

Antibodies Product Datasheet

Product Name:	LIN7B Antibody
Catalog Number:	C49851
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human LIN7B
Concentration:	1 mg/ml
Gene Name:	LIN7B MALS2 VELI2 UNQ3116/PRO10200
Protein Name:	Protein lin-7 homolog B (Lin-7B) (hLin7B) (Mammalian lin-seven protein 2) (MALS-2) (Vertebrate lin-7 homolog 2) (Veli-2) (hVeli2)
Other Name:	
Human Swiss Prot No.:	Q9HAP6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of LIN7B at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year