

Antibodies Product Datasheet

Product Name:	β III Tubulin Antibody
Catalog Number:	C47452
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB: 1:2000-5000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Synthetic Peptide of β III Tubulin
Concentration:	
Gene Name:	TUBB3
Protein Name:	Tubulin beta-3 chain
Other Name:	TUBB3; TUBB4; Tubulin beta-3 chain; Tubulin beta-4 chain; Tubulin beta-III; YM3248
Human Swiss Prot No.:	Q13509
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous β III tubulin protein.
Storage and Stability:	-20°C/1 year