
Antibodies Product Datasheet

Product Name:	ABL1 Antibody
Catalog Number:	R12-2432
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC, ELISA
Dilutions:	ABL1 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ABL1.
Concentration:	1 mg/ml
Gene Name:	ABL1
Other Name:	ABL1; JTK7; Proto-oncogene tyrosine-protein kinase ABL1; c-ABL; p150
Human Swiss Prot No.:	P00519
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	ABL1 Antibodyantibody detects endogenous levels of ABL1.
Storage and Stability:	Stable at -20°C for at least 1 year.