

Antibodies Product Datasheet

Product Name:	CsGalNAcT-2 Antibody
Catalog Number:	C15094
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CSGALNACT2. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	CSGALNACT2
Protein Name:	Chondroitin sulfate N-acetylgalactosaminyltransferase 2
Other Name:	CSGALNACT2; CHGN2; GALNACT2; Chondroitin sulfate N-acetylgalactosaminyltransferase 2; Chondroitin beta-1; 4-N-acetylgalactosaminyltransferase 2; Beta4GalNAcT-2; GalNAcT-2
Human Swiss Prot No.:	Q8N6G5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CsGalNAcT-2 Polyclonal Antibody detects endogenous levels of CsGalNAcT-2 protein.
Storage and Stability:	-20°C/1 year