

Antibodies Product Datasheet

Product Name:	Anti-Human Parathyroid Hormone (ABT215) Mouse Monoclonal Antibody
Catalog Number:	V0431
Clone ID:	ABT215
Reactivity:	Human, Mouse, Rat
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	PTH
Gene Name:	P01270
Other Name:	hPTH, Parathormone, Parathyrin, Parathyroid hormone 1, Parathyroid hormone, Prepro PTH, Preproparathyroid hormone, PTH, PTH1, PTH1 receptor, PTH1R, PTHR, PTHR1, PTHY_HUMAN
Isotype	IgG1, Kappa
Localization:	Parathyroid
Positive Tissue Control:	Cytoplasmic, Membranous
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human Parathyroid Hormone protein.
Storage and Stability:	-20°C/1 year