

Antibodies Product Datasheet

Product Name:	BFSP2 Antibody
Catalog Number:	C49233
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human BFSP2
Concentration:	1 mg/ml
Gene Name:	BFSP2
Protein Name:	Phakinin (49 kDa cytoskeletal protein) (Beaded filament structural protein 2) (Lens fiber cell beaded filament protein CP 47) (CP47) (Lens fiber cell beaded filament protein CP 49) (CP49) (Lens intermediate filament-like light) (LIFL-L)
Other Name:	
Human Swiss Prot No.:	Q13515
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of BFSP2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year