
Antibodies Product Datasheet

Product Name:	hnRNP C1/C2 Antibody
Catalog Number:	R12-2887
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC, ELISA
Dilutions:	hnRNP C1/C2 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human hnRNP C1/C2.
Concentration:	1 mg/ml
Gene Name:	HNRNPC
Other Name:	Heterogeneous nuclear ribonucleoproteins C1/C2; hnRNP C1 / hnRNP C2; HNRNPC; HNRPC
Human Swiss Prot No.:	P07910
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	hnRNP C1/C2 Antibodyantibody detects endogenous levels of hnRNP C1/C2.
Storage and Stability:	Stable at -20°C for at least 1 year.