
Antibodies Product Datasheet

Product Name:	Cortactin (Acetyl-Lys235) Antibody
Catalog Number:	D12055
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized acetyl-peptide derived from the Internal region of human Cortactin around the acetylation site of K235.
Concentration:	1 mg/ml
Gene Name:	CTTN
Protein Name:	Src substrate cortactin
Other Name:	CTTN; EMS1; Src substrate cortactin; Amplixin; Oncogene EMS1
Human Swiss Prot No.:	Q14247
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Acetyl-Cortactin (K235) Polyclonal Antibody detects endogenous levels of Cortactin protein only when acetylation at K235.
Storage and Stability:	-20°C/1 year