

Antibodies Product Datasheet

Product Name:	Tau Antibody
Catalog Number:	C47437
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000 IHC 1:50-300
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic peptide from human protein at AA range: 150-220
Concentration:	1 mg/ml
Gene Name:	MAPT MAPTL MTBT1 TAU
Protein Name:	Microtubule-associated protein tau (Neurofibrillary tangle protein) (Paired helical filament-tau) (PHF-tau)
Other Name:	Microtubule-associated protein tau (Neurofibrillary tangle protein) (Paired helical filament-tau) (PHF-tau)
Human Swiss Prot No.:	P10636
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous Tau
Storage and Stability:	-20°C/1 year