

Antibodies Product Datasheet

Product Name:	EPHA7 (Phospho-Tyr608) Antibody
Catalog Number:	A9517
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human EPHA7 (Phospho Tyr608)
Concentration:	1 mg/ml
Gene Name:	EPHA7 EHK3 HEK11
Protein Name:	EPHA7 (Phospho Tyr608)
Other Name:	Ephrin type-A receptor 7 (EC 2.7.10.1) (EPH homology kinase 3) (EHK-3) (EPH-like kinase 11) (EK11) (hEK11)
Human Swiss Prot No.:	Q15375
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat EPHA7 (Phospho Tyr608)
Storage and Stability:	-20°C/1 year