

Antibodies Product Datasheet

Product Name:	ZN358 Antibody
Catalog Number:	C47299
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human ZN358
Concentration:	1 mg/ml
Gene Name:	ZNF358
Protein Name:	ZN358
Other Name:	Zinc finger protein 358
Human Swiss Prot No.:	Q9NW07
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.160% sodium azide.
Specificity:	This antibody detects endogenous levels of ZN358 at Human
Storage and Stability:	-20°C/1 year