

Antibodies Product Datasheet

Product Name:	IRF9 Antibody
Catalog Number:	C48167
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human IRF9
Concentration:	1 mg/ml
Gene Name:	IRF9 ISGF3G
Protein Name:	Interferon regulatory factor 9 (IRF-9) (IFN-alpha-responsive transcription factor subunit) (ISGF3 p48 subunit) (Interferon-stimulated gene factor 3 gamma) (ISGF-3 gamma) (Transcriptional regulator ISGF3 subunit gamma)
Other Name:	
Human Swiss Prot No.:	Q00978
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of IRF9 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year