

Antibodies Product Datasheet

Product Name:	RIPK2 (Phospho-Ser176) Antibody
Catalog Number:	A8180sec
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized phospho derived from human RIPK2 (Phospho-Ser176)
Concentration:	1 mg/ml
Gene Name:	RIPK2 CARDIAK RICK RIP2 UNQ277/PRO314/PRO34092
Protein Name:	Receptor-interacting serine/threonine-protein kinase 2 (EC 2.7.11.1) (CARD-containing interleukin-1 beta-converting enzyme-associated kinase) (CARD-containing IL-1 beta ICE-kinase) (RIP-like-interacting CLARP kinase) (Receptor-interacting protein 2) (RIP-
Other Name:	Receptor-interacting serine/threonine-protein kinase 2 (EC 2.7.11.1) (CARD-containing interleukin-1 beta-converting enzyme-associated kinase) (CARD-containing IL-1 beta ICE-kinase) (RIP-like-interacting CLARP kinase) (Receptor-interacting protein 2) (RIP-
Human Swiss Prot No.:	O43353
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This detects endogenous levels of RIPK2 (Phospho-Ser176)
Storage and Stability:	-20°C/1 year