

Antibodies Product Datasheet

Product Name:	Ribosomal Protein L10A Antibody
Catalog Number:	C14145
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Ribosomal Protein L10A . at AA range: 30-110
Concentration:	1 mg/ml
Gene Name:	RPL10A
Protein Name:	60S ribosomal protein L10a
Other Name:	RPL10A; NEDD6; 60S ribosomal protein L10a; CSA-19; Neural precursor cell expressed developmentally down-regulated protein 6; NEDD-6
Human Swiss Prot No.:	P62906
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Ribosomal Protein L10A Polyclonal Antibody detects endogenous levels of Ribosomal Protein L10A protein.
Storage and Stability:	-20°C/1 year