

## Antibodies Product Datasheet

<b>Product Name:</b>	ZRANB3 Antibody
<b>Catalog Number:</b>	C48158
<b>Reactivity:</b>	Human, Mouse
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:500-2000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal, Rabbit,IgG
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human ZRANB3
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	ZRANB3
<b>Protein Name:</b>	DNA annealing helicase and endonuclease ZRANB3 (Annealing helicase 2) (AH2) (Zinc finger Ran-binding domain-containing protein 3) [Includes: DNA annealing helicase ZRANB3 (EC 3.6.4.-)Endonuclease ZRANB3 (EC 3.1.-.-)]
<b>Other Name:</b>	
<b>Human Swiss Prot No.:</b>	Q5FWF4
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of ZRANB3 at Human, Mouse
<b>Storage and Stability:</b>	Store at -20° C/1 year