

Antibodies Product Datasheet

Product Name:	KCNA1 Antibody
Catalog Number:	R12-2944
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	KCNA1 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human KCNA1.
Concentration:	1 mg/ml
Gene Name:	KCNA1
Other Name:	Potassium voltage-gated channel subfamily A member 1; Voltage-gated potassium channel subunit Kv1.1; HUK1; HBK1
Human Swiss Prot No.:	Q09470
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	KCNA1 Antibody antibody detects endogenous levels of KCNA1.
Storage and Stability:	Stable at -20°C for at least 1 year.