

Antibodies Product Datasheet

Product Name:	PBK (Phospho-Thr9) Antibody
Catalog Number:	A0710
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PBK/TOPK around the phosphorylation site of Thr9. AA range:1- 50
Concentration:	1 mg/ml
Gene Name:	РВК
Protein Name:	Lymphokine-activated killer T-cell-originated protein kinase
Other Name: Human Swiss Prot No.:	PBK; TOPK; Lymphokine-activated killer T-cell-originated protein kinase; Cancer/testis antigen 84; CT84; MAPKK-like protein kinase; Nori-3; PDZ- binding kinase; Spermatogenesis-related protein kinase; SPK; T-LAK cell-originated protein kinas Q96KB5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-PBK (T9) Polyclonal Antibody detects endogenous levels of PBK protein only when phosphorylated at T9.
Storage and Stability:	-20°C/1 year