
Antibodies Product Datasheet

Product Name:	Claudin 2 mouse mAb(PT0311)
Catalog Number:	N1746
Reactivity:	Human
Applications:	IHC-p,IF(paraffin section),WB
Dilutions:	IHC-p 1:100-300, WB 1:500-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse/IgG2b, Kappa
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human Claudin 2
Concentration:	1 mg/ml
Gene Name:	CLDN2 PSEC0059 SP82 UNQ705/PRO1356
Protein Name:	Claudin 2
Other Name:	Claudin-2 (SP82)
Human Swiss Prot No.:	P57739
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Claudin 2
Storage and Stability:	-20°C/1 year