

Antibodies Product Datasheet

Product Name:	KLF5 Antibody
Catalog Number:	C48053
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KLF5
Concentration:	1 mg/ml
Gene Name:	KLF5 BTEB2 CKLF IKLF
Protein Name:	Krueppel-like factor 5 (Basic transcription element-binding protein 2) (BTE-binding protein 2) (Colon krueppel-like factor) (GC-box-binding protein 2) (Intestinal-enriched krueppel-like factor) (Transcription factor BTEB2)
Other Name:	
Human Swiss Prot No.:	Q13887
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KLF5 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year