

## Antibodies Product Datasheet

|                               |   |
|-------------------------------|---|
| <b>Product Name:</b>          | D55 Antibody  |
| <b>Catalog Number:</b>        | C11501  |
| <b>Reactivity:</b>            | Human, Rat, Mouse   |
| <b>Applications:</b>          | IHC-p, IF(paraffin section), ELISA  |
| <b>Dilutions:</b>             | Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.                            |
| <b>Modification:</b>          | Unmodified/Total  |
| <b>Source:</b>                | Polyclonal Rabbit   |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Immunogen:</b>             | The antiserum was produced against synthesized peptide derived from human TPD52L3. AA range:10-59                     |
| <b>Concentration:</b>         | 1 mg/ml   |
| <b>Gene Name:</b>             | TPD52L3   |
| <b>Protein Name:</b>          | Tumor protein D55   |
| <b>Other Name:</b>            | TPD52L3; Tumor protein D55; hD55; Testis development protein NYD-SP25; Tumor protein D52-like 3                       |
| <b>Human Swiss Prot No.:</b>  | Q96J77  |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.   |
| <b>Specificity:</b>           | D55 Polyclonal Antibody detects endogenous levels of D55 protein.   |
| <b>Storage and Stability:</b> | -20°C/1 year  |