

Antibodies Product Datasheet

Product Name:	E2F-4/5 Antibody
Catalog Number:	C0178
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human E2F4. AA range:51-100
Concentration:	1 mg/ml
Gene Name:	E2F4/E2F5
Protein Name:	Transcription factor E2F4/5
Other Name:	E2F4; Transcription factor E2F4; E2F-4; E2F5; Transcription factor E2F5; E2F-5
Human Swiss Prot No.:	Q16254/Q15329
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	E2F-4/5 Polyclonal Antibody detects endogenous levels of E2F-4/5 protein.
Storage and Stability:	-20°C/1 year